

# Pulmodyne® BiTrac™ NIV

comfortable fit. efficient seal



## Non-Invasive Ventilation

Pulmodyne's BiTrac™ NIV mask is a bilevel compatible and disposable mask with a silicone forehead pad, multi-adjustable fittings using the OmniClip™, and four sizes up to an XL mask. These features provide an efficient seal and comfortable fit for patients.



Complete line  
of masks for  
immediate  
non-invasive  
ventilation

# Ordering information

BiTrac™ NIV Masks product information		Part numbers
<b>BiTrac™ NIV Full Face Mask (Bilevel/CPAP Application)</b> with Anti-Asphyxia Elbow 22mm Female, OmniClip™ with Silicone Forehead Pad and Head Strap		
Adult Extra Large	10 per case	<b>313-9004</b>
Adult Large	10 per case	<b>313-9003</b>
Adult Medium	10 per case	<b>313-9002</b>
Adult Small	10 per case	<b>313-9001</b>
<b>BiTrac™ SE NIV Full Face Mask (Ventilator Application)</b> with Standard Elbow 22mm Female, OmniClip™ with Silicone Forehead Pad and Head Strap		
Adult Extra Large	10 per case	<b>313-9012</b>
Adult Large	10 per case	<b>313-9011</b>
Adult Medium	10 per case	<b>313-9010</b>
Adult Small	10 per case	<b>313-9009</b>
<b>BiTrac™ NIV Nasal Mask (Bilevel/CPAP and Ventilator Applications)</b> with Standard Elbow 22mm Female, Silicone Forehead Pad and Head Strap		
Adult Extra Large	10 per case	<b>313-9016</b>
Adult Large	10 per case	<b>313-9015</b>
Adult Medium	10 per case	<b>313-9014</b>
Adult Small	10 per case	<b>313-9013</b>
Circuits product information		Part numbers
<b>BiTrac™ NIV Circuit (Bilevel/CPAP Applications)</b>		
72" Smooth ID Tubing with swivel connector and exhalation port, cap plug; includes optional 24" tubing with "T" and proximal pressure line	10 per case	<b>313-9074</b>
<b>BiTrac™ NIV Circuit (Bilevel/CPAP Applications)</b>		
72" Smooth ID Tubing with attached filter, swivel connector and exhalation port, cap plug; includes optional 24" tubing with "T" and proximal pressure line	10 per case	<b>313-9075</b>
Accessories		Part numbers
Head Strap with 2 clips (Black)	10 per case	<b>313-9026</b>
Full Cap Head Gear with 2 clips (Black)	10 per case	<b>313-9027</b>
Chin Strap (Black)	10 per case	<b>313-9029</b>

## BiTrac™ NIV Masks



BiTrac™ NIV Full Face Mask



BiTrac™ SE NIV Full Face Mask



BiTrac™ NIV Nasal Mask



BiTrac™ NIV Circuit

For more information please contact us at:

Air Liquide Healthcare, Stratus House, Unit 1 IDA College Park, Blanchardstown Road North, Blanchardstown, Dublin 15, D15 PEC4.

Tel: +353 (0)1 809 1800 Freephone: 1800 24 02 02 (ROI only) 0800 328 5875 (NI only)

Fax: +353 (0)1 829 3966 Email: [healthie@airliquide.ie](mailto:healthie@airliquide.ie) Website: [www.airliquidehealthcare.ie](http://www.airliquidehealthcare.ie)